

DEXCOWIN

DEXCOWIN

IRAYA6

ADX6000FB

DEXCOWIN'S SIGNATURE LINE OF X-RAY SYSTEMS

SPECIFICATIONS

Tube Voltage	40KV~90KV	Exposure	0.05 - 1.5sec
Tube current	1~6 mA	Detector	Flat panel sensor/detector
Focal spot	0.8mm	Power Supply	22.2Vdc
Target angle	16°		
Total Filtration	More than 3.1mmAl(At 90kV) (Inherent Filtration: 1mmAl, Fixed added filter: 1.9mmAl, Collimator: 0.2mmAl) IRAYA6 only : More than 2.5mmAl(At 90kV) (Inherent Filtration: 1mmAl, Fixed added filter: 1.5mmAl)		



FDA CE ISO

DEXCOWIN

IRAY A6

ADX6000FB

DEXCOWIN'S SIGNATURE LINE OF X-RAY SYSTEMS



Medical



Veterinary



Orthopedic



SPECIFICATIONS

Advanced Radiation Safety


- 1/7 lower radiation doses than wall mounted x-ray
- **Internal** Multiple Lead wrapping ensures no radiation leakage

+ **Highlight** Ensures both user & patient safety by reducing radiation exposure

+ **Highlight** Efficient and faster workflow

Battery

- Takes 150 Short at 6L battery
- 6 Cell high capacity battery pack
- Adjustable power



Composition

+ **Standard Package** Main body, 2 Batteries and charger set, Manual CD, Certificate of Conformance, Hard Case

