



BARE LASER FIBER

STRAIGHT
BEAM



Hemorrhoids



Fissure



Dental



Oncology



Fetal
Reduction




ENT



PLDD

Bare Laser Fiber

arphotoniccsfiber.com



Specification

Fiber Core Diameter	200µm	400µm	600µm	800µm	1000µm	
Minimum Bending Radius	10mm	30mm	40mm	60mm	80mm	
Outer Diameter	0.3mm	0.9mm	1.2mm	1.3mm	1.4mm	
Effective Length						2.7m
Wavelength Compatibility						1470, 980 & 1940nm
Maximum Transmission Power						20W
Minimum Transmission Efficiency						96%
Fiber Lifespan						36 months
Sterilization Method						Ethylene Oxide