


Video Rectoscope


DCS-103Pro

Anoscope for examination of internal or External Hemorrhoids, Anal Fissure and Anal Fistula


- Full HD
- Telemedicine system
- Convenient joystick Operation
- Image Conversion
- Image management
- 18W LED Lamp loaded
- Easy Connection with PC/Laptops




Piles





Fissures



Fistula



Pilonidal Sinus



Convenient joystick with each directions leads to freeze, save and zoom images

DCS Pro light source can be easily connected to a computer via USB 2.0. Users can save digital data on a computer and print images



Functions	FULL HD Live Image (16:9) Image Comparison USB2.0 PC Slave Mode Image Adjustment - Brightness, Sharpness, Contrast, Gamma, Exposure - Image Conversion - Selectable color mode(User 1/2/3/Default) - Stored Image replay - System Time setting
Camera Resolution	1920x1080
Video Out	1920x1080p, 60 frame, HDMI
PC Interface	USB 2.0 / Slave Mode Directly connected to the PC application
Internal Storage	32GB
External Storage interface	USB 2.0 Memory Stick /FAT32
Still image Capture	1920x1080 jpg File
User Key Interface	5 Directional Key
Light Source	18W LED lamp, Color temperature (4800. K), Adjustable light intensity
Cable Length	2.0M
Foot switch	for Freeze & Still Capture
AC Input	90~264VAC Max 45Watt Universal AC Input / Full range
Dimension	Width(209mm) x Depth(101mm) x Height(250mm)



Analyticals Medical
Technologies Pvt. Ltd.

AMT, Plot No- A815, MIDC Industrial Area,
Kopar Khairane, Navi Mumbai, Maharashtra 400709

Office : +91-9322673077 | +91-9892703077

Office : +91-8080321777 | +91-9321548914

Web : www.amtmedicals.com

Email : info@amtmedicals.com