



**VIDEO  
RECTOSCOPE**

# ANOCAM PROCTOLOGY

IMAGING CAMERA

**1 ULTRA**  
Portable

**2 AUTO**  
White - Balance

**3 PLUG**  
& Play



Piles



Fissures



Fistula



Pilonidal  
Sinus



**DIAGNOSTIC  
REPORT**



## Features

- Small and Lightweight
- No fiber optic and light source
- Portable
- Works almost with all proctoscope
- Built-in adjustable Illumination
- Integrated software to freeze and store
- USB connectivity
- No power supply or batteries required
- Auto Focus
- Anti-Fog
- Shock Proof
- Image Saving
- Slider to adjust the light illumination
- Capturing buttons (for capturing images)
- USB cable (3 meters)

## Specification

**Name :** Anocam

**Resolution :** 640x720

**Luminescence :** 8 in-built LED lights

**Image trigger control :** Hand-switch

**Focus :** Auto Focus

**White Balance :** Auto white balance

**Frame Rate :** 25FPS

**Viewing Angle :** 62°

**Sensor chip :** 1/6" CMOS

**Video Output :** USB 2.0

**Build :** Polycarbonate enclosure

**Weight :** 400 g

**Kit contents :** Camera, full moon proctoscope, USB key

**Operating System :** Windows 10 and up

Anoscope can be easily connected to a computer via USB 2.0 user can save digital data on a computer & create the report on the go

Anoscope for examination of internal or external hemorroids, anal fissure & anal fistula



The software connects to the camera for image output after installation.



The aluminum box has two layers: top for camera, bottom for a laptop.



All accessories are packed and ready to go

